

Why is the Air/Oxygen Mixer an optional accessory?

The HT50 Ventilator offers the flexibility your patients need in oxygen accessories:

- Oxygen Blending Bag Kit, for Low Pressure (p/n HT600044)
 - Air/Oxygen Entrainment Mixer, for High Pressure (p/n HT460700)
- For homecare patients who do not need more than 21% room air oxygen – a mixer is not required.
 - Patients that need to be ventilated with oxygen but find that a low flow device (such as an oxygen concentrator) is more convenient or economical can use the **Oxygen Blending Bag** to deliver mixed gas.
 - The **Air/Oxygen Entrainment Mixer** is a simple device that attaches to a 50 psi oxygen source (such as an e-cylinder) and provides a high flow mix of 21-100% oxygen by dialing in the appropriate amount.



Oxygen Blending Bag Kit



Air/Oxygen Entrainment Mixer

Breathe easier with the O₂ options available on the HT50.