

TO: Newport Customers

SUBJECT: Silicone Neonatal Test Lung p/n G-118000-0

Newport is pleased to offer a Silicone Neonatal Test Lung p/n G-118000-0 for use in ventilator testing and training.

With a capacity of up to 60 mL and a mean compliance of ~ 0.8 mL/cmH₂O, this pleated, free-style bellows test lung with endotracheal-tube style fitting will allow you to simulate neonatal through small pediatric patients.

Pulling the small tab which projects out from the bottom of the lung allows you to simulate patient triggering on the ventilator by creating negative pressure, similar to that which would be associated with a breathing effort.

To Order

For pricing and to order the Silicone Neonatal Test Lung p/n G-118000-0, please contact Newport Customer Service Dept at +1.714.427.5811 ext 282 or email: customers@ventilators.com.



MKPB012811