

TO: Newport Medical International Sales Managers and Distributors

SUBJECT: e360T Model Has Enhanced Features

The latest model of the e360 Ventilator series, the e360T, includes the addition of an external touch screen display monitor and new software features in addition to features found on the E model.

The external touch screen is a 15" medical grade monitor that allows the user to make changes to the ventilator from the external screen in the same way as using the GUI on the ventilator face panel. This gives the user added convenience and control to use either screen to make changes.

The new software (TSN 3.9) include a PV Maneuver tool, a multiple waveform screen that displays three waveforms and a pop up window that displays I:E Ratio changes in real time.

New Touch Screen !



T Model Software Version TSN 3.9:

- ✓ PV Maneuver - This is a manual slow inflation maneuver that provides an added tool for assessing lung dynamics and setting PEEP. It is simple to perform and can be accessed from the "Technical" screen.
- ✓ Multiple Waveforms on One Screen: Pressure, Volume and Flow waves can be viewed at the same time to make graphics analysis quick and easy.
- ✓ View Real Time Changes to I:E Ratio Settings: While adjusting parameters for Resp Rate, t insp, Flow, Pause or Tidal Volume, a pop-up window shows the I:E Ratio adjustment.
- ✓ All other software features are the same as found in the e360E model software including the "Intelligent Controls System" automated features and Save/Download features.

NOTE: e360 ventilators versions S, P and E cannot be upgraded to the T version. There are hardware changes required for the T version.

Ordering information:

Order part number for the T version: E360T-WW*-IN or E360T-WW*-NA

*Other country specific versions available include CN, ES, PL, PT, RU, SY (symbols).

To order please contact Newport Customer Service at +1.714.427.5811 ext. 282 or email: customers@ventilators.com.

MKPB052610