

To: All Newport Customers

Subject: Discontinued Use of Scroll and Paper in Infant (MR340) & Adult (MR370) Fisher & Paykel Reusable Humidifier Chambers

Over the course of the years, Fisher & Paykel has made technological advancements to their humidifiers. One of the key changes has been to make the heater plates larger, allowing them to output higher levels of humidity. Recently, Fisher & Paykel discovered that their humidifiers work more efficiently when there is no aluminum scroll or paper in the water chamber.

Since the Fisher & Paykel Infant (p/n MR340) & Adult (p/n MR370) Reusable Humidifier Chambers provide more optimal humidity when used without scrolls and paper and therefore it is no longer recommended to use the aluminum scroll and paper, Fisher & Paykel has stopped including these items in the packaging of the Infant and Adult Reusable Humidifier Chambers.



The Fisher & Paykel Infant (p/n MR340) & Adult (p/n MR370) Reusable Humidifier Chambers that are shipped from Newport no longer contain aluminum scrolls or absorbent paper. Newport has also stopped ordering the replacement scrolls and paper.

We have a limited quantity of Absorbent Paper Refills available. When these supplies are gone, we will discontinue selling them.

P/N 900MR066 Absorbent Paper Refill Pack for MR370

P/N 900MR065 Absorbent Paper Refill Pack for MR340

For more information:

Please contact Newport Customer Service at +1.714.427.5811 ext. 282 or email: customers@ventilators.com.

MKPB110911