



November 30, 2010

Re: HT70 Use During Travel by Air

To Whom It May Concern:

This is to confirm that the HT70 series of ventilators utilizes a sealed battery that complies with non-spillable requirements such as 49 CFR part 173.159. The HT70 has also been tested to and passed DO-160D, section 21 specifically concerning emissions. The HT70 is cleared for use on commercial air flights during take off, landing and flight.

The HT70 has been used on commercial air flights all over the world to transport people needing ventilator support. It is also designed for use in helicopter and fixed wing medical life flight systems.

In view of the design and documented verification that the HT70 meets these requirements, the HT70 series of ventilators are appropriate to use for travel on commercial airlines.

Sincerely,

A handwritten signature in blue ink that reads "Jim Seeley".

Jim Seeley, RRT
Product Manager
Newport Medical Instruments
714-427-5811, ext. 149
email: JSeeley@NewportMedical.com

Newport Medical Instruments, Inc.

1620 Sunflower Avenue, Costa Mesa, CA 92626 USA • +1.714.427.5811 • fax: +1.714.427.0489 •
www.ventilators.com